

Finance Options

LactoScope Filter



The LactoScope Filter analyser is available on a finance scheme^a allowing you to avoid capital outlay and spread your payments over 5 years. This provides you with a brand new LactoScope with minimal up-front expenditure and spreads the cost over time via fixed payments. The scheme will finance everything related to the purchase and installation of the analyser and includes annual preventative maintenance visits.

Payments are 100% deductible against taxable profits and can be structured to fit your requirements. The scheme reduces the drain on your working capital and is generally not classified as debt but as a running cost or expense.

A New Analyser Every 5 Years

After 5 years you have the option to buy the analyser for a nominal fee or renew the scheme and receive a brand new replacement analyser.

LactoScope Filter C3+		
Term	Monthly Payment	Monthly Payment After Tax Relief
60 Months	£559.17 ^b	£447.34 ^c

LactoScope Filter

The Delta Instruments LactoScope Filter is a stand-alone versatile infrared dairy analyser ideal for liquid dairy and cheese processors who require simple, accurate and cost effective routine sample analysis.

- Up to 150 samples per hour
- Automatic sample heating, zeroing and cleaning with 5 litre zeroing and cleaning tanks
- Generic cleaning/zero solutions and no standardisation solution required
- 20 product menu
- Low maintenance costs
- Freezing point depression (option)

LactoScope Filter Specification	
Fat	0-60%
Protein	0-15%
Lactose	0-20%
Total Solids	0-60%
Solids-Non-Fat	0-20%
FPD	450°C - 550°C M (optional)
Accuracy	Typically <1% CV. FPD 0.005°C
Repeatability	Typically <0.25% CV
Sample Volume	Typically 8 ml
Sample Temperature	2-42°C

^a Subject to credit rating

^b Not including VAT

^c Assumes corporation tax at 20%

Prices shown in this document are correct as of 23/03/2017 and are subject to change.