

On-Site LactoScope FTIR Trials



DeltaInstruments
An Advanced Instruments Company

LactoScope FTIR

**Measure Fat, Protein, Total Solids,
Solids non Fat, Casein and FPD in:**

- ❖ **Milk**
- ❖ **Cream**
- ❖ **Cheese**
- ❖ **Whey**
- ❖ **Yogurt**
- ❖ **Ice Cream**

The LactoScope FTIR is a robust and easy to use analyser providing exceptional accuracy and repeatability. Measuring up to 120 samples per hour it is suitable for highly viscous products and features automated zeroing, cleaning and sample heating.

Our demonstration units are available for on-site trials, giving you the opportunity to compare it side by side with your current analyser. If you would like to proceed with a trial we will:

- ❖ Deliver and install the LactoScope
- ❖ Calibrate the analyser using DRRR Standardon standards
- ❖ Provide full training for your laboratory staff
- ❖ Leave the LactoScope at your site for up to 4 weeks
- ❖ Review performance and findings at the end of the trial period
- ❖ Uninstall and remove the analyser

To arrange a trial please call 01342 820820 or e-mail sales@qclscientific.com.