

SpectraAnalyzer

SAMPLE PRESENTATION ACCESSORIES





1

Closed Cup

for the analysis of powders:
e.g. milk powder, starch,
ground coffee and flour.



2

Open Cup

for the analysis of
pastes: e.g. butter
and cream cheese.



5

Slide Cup

for the analysis of products
with low water content:
e.g. oil and glycerine.



6

Italian Cup

for the analysis of chocolate mass
and for clinical analysis
of human faeces.



9

Liquid Cell Assy Top Window

for the analysis of liquid samples
for the SpectraAnalyzer Top Window
only: e.g. oil, glycerine, marinades
and ketchup.



3

Open Cup, Deep

for the analysis of olive paste,
olive cake and minced meat.



4

Viscous Cup

for the analysis of viscous samples:
e.g. honey, ketchup and
condensed milk.



7

Cheese Cup

for the analysis of
hard cheese.



8

Glass Cup TW

for the analysis of powders for the
SpectraAlyzer Top Window only:
e.g. milk powder, starch, ground
coffee and flour.



Cups 1 to 7 are sample presentation cups for the SpectraAlyzer® drawer models,
cups 8 and 9 are for the use with the SpectraAlyzer® Top Window model only.



Product name and item number

Product name	Item no.
1. Closed Cup	10206
2. Open Cup	10201
3. Open Cup, Deep	10202
4. Viscous Cup	10210
5. Slide Cup	10203
6. Italian Cup	10209
7. Cheese Cup Assy	10260
8. Glass Cup, TW	10208
9. Liquid Cell Assy Top Window	10265

